


Recomended Basic Laryngology Set

- LN-3503 Laryngoscope Rod and Chest Support,**
- LN-3507 Laryngoscope plate small** ,18 cm, for children, it can be used with LN-3502.
- LN-3504 Laryngoscope plate small** ,18 cm, for adults, it can be used with LN-3502.
- LN-3505 Laryngoscope plate medium** ,18 cm, for adults, it can be used with LN-3502.
- LN-3506 Laryngoscope plate large** ,18 cm, for adults, it can be used with LN-3502.
- LN-3508 Ferhan Öz type Laryngoscope plate**
length 15 cm, for adults,distal 18 mm, it can be used with LN-3502.
- LN-3510 Micro-laryngeal cup shaped Forceps**, straight, Ø 2 mm, 25 cm long.
- LN-3513 Micro-laryngeal cup shaped Forceps**, angled left, Ø 2 mm, 25 cm long.
- LN-3520 Micro-laryngeal Scissors**, straight, 25 cm long.
- LN-3522 Micro-laryngeal Scissors**, angled left, 25 cm long.
- LN-3525 Micro-laryngeal Scissors**, angled right, 25 cm long.
- LN-3527 Micro-laryngeal Scissors**, angled up, 25 cm long.





LN-3559 Micro-laryngeal full set,
7 instruments with storage
case and handle.

**LN-3515 Micro-laryngeal cup
shaped Forceps,** angled
right, Ø 2 mm, 25 cm long.

**LN-3517 Micro-laryngeal cup
shaped Forceps,** up,
Ø 2 mm, 25 cm long.

**LN-3530 Micro-laryngeal
Grasping Forceps,**
Straight, 25 cm long.

**LN-3533 Micro-laryngeal
Grasping Forceps,**
angled left, 25 cm long.

**LN-3535 Micro-laryngeal
Grasping Forceps,**
angled right, 25 cm long.

**LN-3540 Micro-laryngeal cup
shaped (Punch)
Forceps,** for biopsy universal,
25 cm long.

**LN-3543&LN-3545 Suction Tubes
with finger cut-off,**
25 cm long.

**LN-3538 Ferhan Öz type basket
Forceps,** left up tips,
23 cm working length .

**LN-3539 Ferhan Öz type basket
forceps,** right up tips,
23 cm working length .

LN-3518 Fiber Optik Light Carrier
for use with Laryngoscope,
135 mm length. it can be
used with LN-3504.LN-3505.
LN-3506.LN-3508. LN-3509.

LN-3519 Fiber Optik Light Carrier,
for use with Laryngoscope,
85 mm length. it can be
used with LN-3504.LN-3505.
LN-3506.LN-3508. LN-3509.

